

# BUKETS FOR PUNCTURE-CUTTING WASTE

steriale® | BOX

## CHARACTERISTICS



Volume:  
Material:  
Use:  
Sterilized method:  
Characteristics:

Dimensions:

Net weight:

1, 2.5, 5, 10 Liters  
paperboard/ plastic cap  
puncture-cutting waste  
By steam pressure water/AUTOCLAVE  
ONLY ONE USE SUPPLY It is used together  
with reusable containers Steriale Containers

1L Ø100mm x 120mm / 2.5L Ø200mm x 160mm  
5L Ø230mm x 240mm / 10L Ø230mm x 360mm

1L (Pack of 50 U.) 100 Gr /  
2.5L (Pack of 50 U.) 200 Gr /  
3L (Pack of 50 U.) 300 Gr /  
5L (Pack of 50 U.) 400 Gr /  
10L (Pack of 50 U.) 700 Gr /



00066 SP1

barcode on cartridges for traceability, change color  
with sterilization.

# TRIPLE LAYER BAGS FOR SOLID WASTE

## CHARACTERISTICS



Volume:	16, 36 Liters.
Material:	Low density polypropylene plastic, paper, heat shrinkable. No puncture-cutting waste.
Use:	By steam pressure water/AUTOCLAVE.
Sterilized method:	Only one use supply.
Characteristics:	It is used together with reusable containers. Steriale Containers- Hermetics.
Dimensions:	16L Ø220mm x 400mm / 36L Ø280mm x 550mm /
Net weight:	16L (Packs of 100 U.) 2 Kg / 36L (Packs of 100 U.) 4 Kg /
Colours:	Red, Brown or White.



barcode on cartridges for traceability, change color with sterilization.