

PATENTED LID
OPENING SYSTEM
AT THE END OF THE **STERILIZATION**

steriale | **CONTAINERS**



CHARACTERISTICS

Volume:	1, 2.5, 5, 10 Litros
Material:	stainless steel and antipyretic plastic
Use:	By steam/AUTOCLAVE
Sterilized method:	By steam pressure water/ AUTOCLAVE
Characteristics:	REUSABLE it is used together with disposable buckets
Dimensions:	Steriale Box Ø 100mm x 120mm / Ø 100mm x 240mm / Ø 100mm x 350mm / Ø 200mm x 160mm/ Ø 200mm x 320mm /
Net weight:	1L 100 Gr / 2L 200 Gr / 3L 300 Gr / 5L 400 Gr/ 10L 700 Gr/
Colours:	Red y Yellow

PATENTED LID
OPENING SYSTEM
AT THE END OF THE **STERILIZATION**

CHARACTERISTICS



Volume:
Material:
Use:
Sterilized method:
Characteristics:

Dimensions:

Net weight:

Colours:

16, 36 Liters
Stainless steel and antipyretic plastic
By steam/ AUTOCLAVE
By steam pressure water/ AUTOCLAVE
REUSABLE.

It is used together with disposable bags.
STERIALE BAGS.

360mm x Ø 225mm /
530mm x Ø 300mm /

16L 1.6 Kg / 36L 3 Kg /

Red, Black and Green.